

Product name

SIK1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SIK1 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human colon cancer tissue, human stomach tissue.

Preparation

This antibody was obtained by immunization of SIK1 recombinant protein (Accession Number: NM_173354). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC


Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations

Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:115229(SIK1 Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:115229(SIK1 Antibody) at dilution of 1:50 (under 10x lens)