

# SIK1 antibody

Catalog Number: 115229

## Product name

SIK1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SIK1 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human colon cancer tissue, human stomach tissue.

## Preparation

This antibody was obtained by immunization of SIK1 recombinant protein (Accession Number: NM\_173354). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:115229(SIK1 Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:115229(SIK1 Antibody) at dilution of 1:50 (under 10x lens)