

SIRP BETA 1 antibody



Catalog Number: 115239

Product name

SIRP BETA 1 antibody

Specificity

Human; other species not tested.

Antibody description

SIRP BETA 1 Rabbit Polyclonal antibody. Positive WB detected in U-937 cells, A375 cells, PC-3 cells. Positive IHC detected in human lung cancer tissue. Observed molecular weight by Western-blot: 43 kDa

Preparation

This antibody was obtained by immunization of SIRP BETA 1 recombinant protein (Accession Number: NM_001083910). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:115239(SIRPB1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:115239(SIRPB1 antibody) at dilution of 1:50



(under 40x lens)

U-937 cells were subjected to SDS PAGE followed by western blot with Catalog No:115239(SIRPB1 Antibody) at dilution of 1:1000