

SIRT4 antibody



Catalog Number: 115245

Product name

SIRT4 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

SIRT4 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, HeLa cells, NIH/3T3 cells. Positive IHC detected in human skeletal muscle tissue. Observed molecular weight by Western-blot: ~41 kDa

Preparation

This antibody was obtained by immunization of SIRT4 recombinant protein (Accession Number: XM_006719309). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:

WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:115245(SIRT4 Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue slide using Catalog No:115245(SIRT4 Antibody) at dilution of 1:50 (under 40x lens)