

SKIL antibody

Catalog Number: 115257

Product name

SKIL antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SKIL Rabbit Polyclonal antibody. Positive WB detected in A431 cells, HepG2 cells, human skeletal muscle tissue. Positive IP detected in HepG2 cells. Observed molecular weight by Western-blot: 77 kDa

Preparation

This antibody was obtained by immunization of SKIL recombinant protein (Accession Number: NM_005414). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:115257(SKIL antibody) at dilution of 1:300



IP Result of anti-SKIL (IP:Catalog No:115257, 4ug; Detection:Catalog No:115257 1:500) with HepG2 cells lysate 1600ug.