

# SHC4 antibody



Catalog Number: 115268

## Product name

SHC4 antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

SHC4 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, HepG2 cells, human brain tissue. Observed molecular weight by Western-blot: 70-75kd

## Preparation

This antibody was obtained by immunization of SHC4 recombinant protein (Accession Number: NM\_203349). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:115268(SHC4 antibody) at dilution of 1:300