

SHC4 antibody



Catalog Number: 115268

Product name

SHC4 antibody

Specificity

human, mouse; other species not tested.

Antibody description

SHC4 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, HepG2 cells, human brain tissue. Observed molecular weight by Western-blot: 70-75kd

Preparation

This antibody was obtained by immunization of SHC4 recombinant protein (Accession Number: NM_203349). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:115268(SHC4 antibody) at dilution of 1:300