

# SHISA4 antibody

Catalog Number: 115274

## Product name

SHISA4 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

SHISA4 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: 25-27 kDa

## Preparation

This antibody was obtained by immunization of SHISA4 recombinant protein (Accession Number: XM\_017000419). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115274(SHISA4 antibody) at dilution of 1:2000



Immunofluorescent analysis of HEK-293 cells using Catalog No:115274(SHISA4 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG