

# SLC13A4 antibody



Catalog Number: 115303

## Product name

SLC13A4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SLC13A4 Rabbit Polyclonal antibody. Positive IF detected in HEK-293 cells. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 69 kDa

## Preparation

This antibody was obtained by immunization of SLC13A4 recombinant protein (Accession Number: NM\_012450). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115303(SLC13A4 antibody) at dilution of 1:500



Immunofluorescent analysis of (10% Formaldehyde) fixed HEK-293 cells using Catalog No:115303(SLC13A4 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)