

# SLC16A12 antibody



Catalog Number: 115304

## Product name

SLC16A12 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SLC16A12 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue. Positive WB detected in mouse kidney tissue. Observed molecular weight by Western-blot: 53kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_213606).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol  
pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations




mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:115304(SLC16A12 antibody) at dilution of



1:100

Immunohistochemical of paraffin-embedded human testis using Catalog No:115304(SLC16A12 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human testis using Catalog No:115304(SLC16A12 antibody) at dilution of 1:50 (under 40x lens)