

SLC24A6 antibody



Catalog Number: 115319

Product name

Polyclonal

SLC24A6 antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

SLC24A6 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, human kidney tissue. Observed molecular weight by Western-blot: 64-66 kDa

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of SLC24A6 recombinant protein (Accession Number: NM_024959). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:115319(SLC24A6 antibody) at dilution of 1:300