


SLC25A33 antibody



Catalog Number: 115333

Product name	Polyclonal
SLC25A33 antibody	Ig Type
Specificity	Rabbit IgG
Human; other species not tested.	Applications
Antibody description	ELISA, WB
SLC25A33 Rabbit Polyclonal antibody. Positive WB detected in U2Os cells. Observed molecular weight by Western-blot: 35-40 kDa	Dilutions
Preparation	Recommended Dilution:
This antibody was obtained by immunization of SLC25A33 recombinant protein (Accession Number: NM_032315). Purification method: Antigen affinity purified.	WB: 1:200-1:2000
Formulation	Validations
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.	
Storage	WB result of anti-SLC25A33 (Catalog No:115333) by Dr. Geraldine Kelly, University College Cork, Ireland.
Store at -20°C. DO NOT ALIQUOT	
Clonality	