

SLC25A34 antibody

Catalog Number: 115334

Product name

SLC25A34 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SLC25A34 Rabbit Polyclonal antibody. Positive IP detected in HepG2 cells. Positive WB detected in HepG2 cells, human liver tissue. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of SLC25A34 recombinant protein (Accession Number: NM_207348). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:115334(SLC25A34 antibody) at dilution of



1:500

IP Result of anti-SLC25A34 (IP:Catalog No:115334, 4ug; Detection:Catalog No:115334 1:500) with HepG2 cells lysate 1600ug.