

SLC25A45 antibody



Catalog Number: 115336

Product name

SLC25A45 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SLC25A45 Rabbit Polyclonal antibody. Positive IF detected in SH-SY5Y cells. Positive WB detected in mouse skeletal muscle tissue, SH-SY5Y cells. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of SLC25A45 recombinant protein (Accession Number: BC041100). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:115336(SLC25A45 antibody) at dilution of



Immunofluorescent analysis of SH-SY5Y cells using Catalog No:115336(SLC25A45 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)