

SLC5A8 antibody



Catalog Number: 115346

Product name

SLC5A8 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SLC5A8 Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive WB detected in rat kidney tissue, HEK-293 cells, human brain tissue, mouse kidney tissue. Observed molecular weight by Western-blot: 67kd

Preparation

This antibody was obtained by immunization of SLC5A8 recombinant protein (Accession Number: NM_145913). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



Immunofluorescent analysis of Hela cells, using SLC5A8 antibody Catalog No:115346 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

rat kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:115346(SLC5A8 Antibody) at dilution of 1:1000