

# SLC6A8 antibody



Catalog Number: 115347

## Product name

SLC6A8 antibody

## Specificity

Human; other species not tested.

## Antibody description

SLC6A8 Rabbit Polyclonal antibody. Positive WB detected in HUVEC cells. Observed molecular weight by Western-blot: 55-65 kDa

## Preparation

This antibody was obtained by immunization of SLC6A8 recombinant protein (Accession Number: NM\_005629). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HUVEC cells were subjected to SDS PAGE followed by western blot with Catalog No:115347(SLC6A8 Antibody) at dilution of 1:600