

SLC7A4 antibody



Catalog Number: 115348

Product name

SLC7A4 antibody

Specificity

Human; other species not tested.

Antibody description

SLC7A4 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue. Positive IF detected in Hela cells. Observed molecular weight by Western-blot: 68kd

Preparation

This antibody was obtained by immunization of SLC7A4 recombinant protein (Accession Number: NM_004173). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:115348(SLC7A4 antibody) at dilution of 1:500



Immunofluorescent analysis of Hela cells, using SLC7A4 antibody Catalog No:115348 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).