

SLC9A11 antibody



Catalog Number: 115350

Product name

SLC9A11 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

SLC9A11 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HepG2 cells. Observed molecular weight by Western-blot: 135-140kd

Preparation

This antibody was obtained by immunization of SLC9A11 recombinant protein (Accession Number: NM_178527). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:115350(SLC9A11 antibody) at dilution of 1:500