

SLFN1 antibody



Catalog Number: 115356

Product name

SLFN1 antibody

WB: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human; other species not tested.

Validations

Antibody description


SLFN1 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells. Positive IHC detected in human kidney tissue, human brain tissue, human liver tissue, human ovary tissue, human placenta tissue, human skin tissue, human spleen tissue, human testis tissue. Observed molecular weight by Western-blot: 46 kDa



Preparation

This antibody was obtained by immunization of SLFN1 recombinant protein (Accession Number: XM_006710432). Purification method: Antigen affinity purified.

COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:115356(SLFN1 antibody) at dilution of 1:800



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Immunohistochemical of paraffin-embedded human kidney using Catalog No:115356(SLFN1 antibody) at dilution of 1:100 (under 10x lens)

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:



Immunohistochemical of paraffin-embedded human kidney using Catalog No:115356(SLFN1

SLFNL1 antibody



Catalog Number: 115356

antibody) at dilution of 1:100 (under 40x lens)