

# SLK antibody

Catalog Number: 115361

## Product name

SLK antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

SLK Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, human brain tissue, MCF7 cells, mouse liver tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 143-160kd,210-220kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_014720). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115361(SLK antibody) at dilution of 1:2000



Immunofluorescent analysis of HepG2 cells, using SLK antibody Catalog No:115361 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).