

SLC36A4 antibody



Catalog Number: 115369

Product name

SLC36A4 antibody

Specificity

human; other species not tested.

Antibody description

SLC36A4 Rabbit Polyclonal antibody. Positive WB detected in A549 cells, HEK-293 cells, PC-3 cells.

Preparation

This antibody was obtained by immunization of SLC36A4 recombinant protein (Accession Number: NM_152313). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:115369(SLC36A4 antibody) at dilution of 1:1000