

# SLC36A4 antibody



Catalog Number: 115369

## Product name

SLC36A4 antibody

## Specificity

human; other species not tested.

## Antibody description

SLC36A4 Rabbit Polyclonal antibody. Positive WB detected in A549 cells, HEK-293 cells, PC-3 cells.

## Preparation

This antibody was obtained by immunization of SLC36A4 recombinant protein (Accession Number: NM\_152313). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:115369(SLC36A4 antibody) at dilution of 1:1000