

SLC45A2 antibody

Catalog Number: 115379

Product name

SLC45A2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SLC45A2 Rabbit Polyclonal antibody. Positive IP detected in HL-60 cells. Positive WB detected in HL-60 cells, MCF7 cells. Observed molecular weight by Western-blot: 50-55 kDa

Preparation

This antibody was obtained by immunization of SLC45A2 recombinant protein (Accession Number: NM_001297417). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:115379(SLC45A2 antibody) at dilution of 1:1000



IP Result of anti-SLC45A2 (IP:Catalog No:115379, 4ug; Detection:Catalog No:115379 1:500) with HL-60 cells lysate 1600ug.