

# SMCHD1 antibody



Catalog Number: 115384

## Product name

SMCHD1 antibody

## Specificity

Human; other species not tested.

## Antibody description

SMCHD1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells. Observed molecular weight by Western-blot: 220 kDa

## Preparation

This antibody was obtained by immunization of SMCHD1 recombinant protein (Accession Number: BC035774). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115384(SMCHD1 Antibody) at dilution of 1:600