

SMAGP antibody



Catalog Number: 115424

Product name

SMAGP antibody

Specificity

Human; other species not tested.

Antibody description

SMAGP Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells. Observed molecular weight by Western-blot: 25-28kd

Preparation

This antibody was obtained by immunization of SMAGP recombinant protein (Accession Number: BC003379). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:115424(Loc57228 antibody) at dilution of 1:600