

# SNRNP200 antibody

Catalog Number: 115454

## Product name

SNRNP200 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

SNRNP200 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, Jurkat cells, NIH/3T3 cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 200-245 kDa

## Preparation

This antibody was obtained by immunization of SNRNP200 recombinant protein (Accession Number: BC001417). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115454(SNRNP200 antibody) at dilution of



1:1000

IP Result of anti-SNRNP200 (IP:Catalog No:115454, 4ug; Detection:Catalog No:115454 1:1000) with HeLa cells lysate 4000ug.