

# SNRPD2 antibody



Catalog Number: 115461

## Product name

SNRPD2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SNRPD2 Rabbit Polyclonal antibody. Positive IHC detected in human skin tissue, human lung cancer tissue. Positive WB detected in A2780 cells, A549 cells, HL-60 cells. Observed molecular weight by Western-blot: 16.5 kDa

## Preparation

This antibody was obtained by immunization of SNRPD2 recombinant protein (Accession Number: NM\_001384647). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



A2780 cells were subjected to SDS PAGE followed by western blot with Catalog No:115461(SNRPD2 antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human skin tissue slide using Catalog No:115461(SNRPD2 Antibody) at dilution of 1:50