

# SNX11 antibody



Catalog Number: 115470

## Product name

SNX11 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SNX11 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of SNX11 recombinant protein (Accession Number: NM\_152244). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115470(SNX11 Antibody) at dilution of 1:500