

# SNX15 antibody



Catalog Number: 115471

---

## Product name

SNX15 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

SNX15 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, COLO 320 cells, human testis tissue. Observed molecular weight by Western-blot: 50kd and 38kd

## Preparation

This antibody was obtained by immunization of SNX15 recombinant protein (Accession Number: NM\_013306). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:115471(SNX15 antibody) at dilution of 1:400