

# SNX19 antibody

Catalog Number: 115474

## Product name

SNX19 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SNX19 Rabbit Polyclonal antibody. Positive IP detected in Raji cells. Positive WB detected in Raji cells, K-562 cells. Observed molecular weight by Western-blot: 109-110 kDa

## Preparation

This antibody was obtained by immunization of SNX19 recombinant protein (Accession Number: NM\_001347921). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP

## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:115474(SNX19 antibody) at dilution of 1:500



IP Result of anti-SNX19 (IP:Catalog No:115474, 4ug; Detection:Catalog No:115474 1:500) with Raji cells lysate 1600ug.