

SNX21 antibody



Catalog Number: 115476

Product name

SNX21 antibody

Specificity

Human; other species not tested.

Antibody description

SNX21 Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue.

Preparation

This antibody was obtained by immunization of SNX21 recombinant protein (Accession Number: NM_152897). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:115476(SNX21 Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:115476(SNX21 Antibody) at dilution of 1:200 (under 40x lens).