

SNX3 antibody



Catalog Number: 115480

Product name

SNX3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SNX3 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells, A549 cells, HeLa cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 16-20 kDa

Preparation

This antibody was obtained by immunization of SNX3 recombinant protein (Accession Number: NM_003795). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:115480(SNX3 antibody) at dilution of 1:500



Immunofluorescent analysis of HeLa cells using Catalog No:115480(SNX3 Antibody) at dilution of 1:25 and Alexa Fluor 594-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)