

SNX4 antibody



Catalog Number: 115482

Product name

SNX4 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

SNX4 Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 50kd

Preparation

This antibody was obtained by immunization of SNX4 recombinant protein (Accession Number: NM_003794). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:115482(SNX4 antibody) at dilution of 1:300