

# SOCS5 antibody



Catalog Number: 115490

## Product name

SOCS5 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SOCS5 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells, HepG2 cells. Positive IF detected in BxPC-3 cells. Observed molecular weight by Western-blot: 61 kDa

## Preparation

This antibody was obtained by immunization of S<sup>o</sup>CS5 recombinant protein (Accession Number: NM\_144949). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:115490(SocS5 antibody) at dilution of 1:200



Immunofluorescent analysis of BxPC-3 cells using Catalog No:115490(SocS5 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)