

# SOD3 antibody



Catalog Number: 115494

## Product name

SOD3 antibody

## Specificity

Human; other species not tested.

## Antibody description

SOD3 Rabbit Polyclonal antibody. Positive IP detected in human placenta tissue. Positive WB detected in human heart tissue, human placenta tissue. Observed molecular weight by Western-blot: 25-33 kDa

## Preparation

This antibody was obtained by immunization of SOD3 recombinant protein (Accession Number: NM\_003102). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:500-1:5000

## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:115494(SOD3 antibody) at dilution of 1:1000



IP Result of anti-SOD3 (IP:Catalog No:115494, 4ug; Detection:Catalog No:115494 1:1000) with human placenta tissue lysate 2800ug.