

SOD3 antibody

Catalog Number: 115494

Product name

SOD3 antibody

Specificity

Human; other species not tested.

Antibody description

SOD3 Rabbit Polyclonal antibody. Positive IP detected in human placenta tissue. Positive WB detected in human heart tissue, human placenta tissue. Observed molecular weight by Western-blot: 25-33 kDa

Preparation

This antibody was obtained by immunization of SOD3 recombinant protein (Accession Number: NM_003102). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:500-1:5000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:115494(SOD3 antibody) at dilution of 1:1000



IP Result of anti-SOD3 (IP:Catalog No:115494, 4ug; Detection:Catalog No:115494 1:1000) with human placenta tissue lysate 2800ug.