

# SOHLH1 antibody



Catalog Number: 115497

## Product name

SOHLH1 antibody

## Specificity

Human; other species not tested.

## Antibody description

SOHLH1 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue.

## Preparation

This antibody was obtained by immunization of SOHLH1 recombinant protein (Accession Number: NM\_001012415). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:115497(SOHLH1 Antibody) at dilution of 1:100 (under 10x lens).



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:115497(SOHLH1 Antibody) at dilution of 1:100 (under 40x lens).