

SORBS2 antibody



Catalog Number: 115502

Product name

SORBS2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

SORBS2 Rabbit Polyclonal antibody. Positive IP detected in mouse heart tissue. Observed molecular weight by Western-blot: 70-75 kDa

Preparation

This antibody was obtained by immunization of SORBS2 recombinant protein (Accession Number: XM_017008771). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:200-1:2000

Validations



IP Result of anti-SORBS2 (IP:Catalog No:115502, 4ug; Detection:Catalog No:115502 1:500) with mouse heart tissue lysate 4000ug.