

# SOX12 antibody



Catalog Number: 115510

## Product name

Polyclonal

SOX12 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse; other species not tested.

## Applications

## Antibody description

ELISA, WB

SOX12 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, mouse brain tissue. Observed molecular weight by Western-blot: 35-40 kDa

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of SOX12 recombinant protein (Accession Number: NM\_006943). Purification method: Antigen Affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115510(SOX12 Antibody) at dilution of 1:1000