

SP140L antibody

Catalog Number: 115524

Product name

SP140L antibody

Specificity

Human; other species not tested.

Antibody description

SP140L Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 67-77 kDa

Preparation

This antibody was obtained by immunization of SP140L recombinant protein (Accession Number: BC004921). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115524(SP140L Antibody) at dilution of 1:600



IP Result of anti-SP140L (IP:Catalog No:115524, 4ug; Detection:Catalog No:115524 1:500) with HeLa cells lysate 3200ug.