

SPACA3 antibody



Catalog Number: 115528

Product name

SPACA3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SPACA3 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, mouse testis tissue. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: 28-32 kDa, 46-50 kDa

Preparation

This antibody was obtained by immunization of SPACA3 recombinant protein (Accession Number: NM_001317225). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:115528(SPACA3 antibody) at dilution of 1:300



Immunofluorescent analysis of HEK-293 cells using Catalog No:115528(SPACA3 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)