

# SPAG11A antibody



Catalog Number: 115529

## Product name

SPAG11A antibody

WB: 1:500-1:5000

IHC: 1:20-1:200


## Specificity

Human; other species not tested.

## Validations

## Antibody description

SPAG11A Rabbit Polyclonal antibody. Positive WB detected in human testis tissue, human brain tissue. Positive IHC detected in human testis tissue, human kidney tissue, human lung tissue, human ovary tissue, human placenta tissue, human skin tissue, human spleen tissue. Observed molecular weight by Western-blot: 20kd




## Preparation

This antibody was obtained by immunization of SPAG11A recombinant protein (Accession Number: BC058833). Purification method: Antigen affinity purified.

human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115529(SPAG11A antibody) at dilution of 1:300

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT


## Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human testis using Catalog No:115529(SPAG11A antibody) at dilution of 1:50 (under 10x lens)

## Ig Type

Rabbit IgG



## Applications

ELISA, WB, IHC

## Dilutions

Recommended Dilution:

Immunohistochemical of paraffin-embedded

# SPAG11A antibody



Catalog Number: 115529

---

human testis using Catalog No:115529(SPAG11A antibody) at dilution of 1:50 (under 40x lens)