

SPATA16 antibody



Catalog Number: 115549

Product name

SPATA16 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SPATA16 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue.

Preparation

This antibody was obtained by immunization of SPATA16 recombinant protein (Accession Number: NM_031955). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:115549 (SPATA16 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:115549 (SPATA16 Antibody) at dilution of 1:50 (under 40x lens)