

# SPHK1 antibody

Catalog Number: 115565

## Product name

SPHK1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

SPHK1 Rabbit Polyclonal antibody. Positive WB detected in L02 cells. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 42 kDa

## Preparation

This antibody was obtained by immunization of SPHK1 recombinant protein (Accession Number: NM\_001355139). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:115565(SPHK1 antibody) at dilution of 1:100



(under 10x lens)

L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:115565(SPHK1 Antibody) at dilution of 1:300