

SPIN2B antibody



Catalog Number: 115572

Product name

SPIN2B antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SPIN2B Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 29 kDa

Preparation

This antibody was obtained by immunization of SPIN2B recombinant protein (Accession Number: XM_017029552). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115572 (SPIN2B antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:115572 (SPIN2B antibody) at dilution of 1:50 (under 25x lens)