

Product name

SRPR antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SRPR Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 65-70kd

Preparation

This antibody was obtained by immunization of SRPR recombinant protein (Accession Number: NM_003139). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:115597(SRPR antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:115597(SRPR antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:115597(SRPR antibody) at dilution of 1:50 (under 40x lens)