

# SSTR5 antibody

Catalog Number: 115615

---

## Product name

SSTR5 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

SSTR5 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue. Observed molecular weight by Western-blot: 50 kDa, 65 kDa

## Preparation

This antibody was obtained by immunization of SSTR5 recombinant protein (Accession Number: NM\_001053). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:115615(SSTR5 antibody) at dilution of 1:1200