

SSX1 antibody



Catalog Number: 115617

Product name

SSX1 antibody

Specificity

Human; other species not tested.

Antibody description

SSX1 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 21-25 kDa

Preparation

This antibody was obtained by immunization of SSX1 recombinant protein (Accession Number: NM_001278691). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:115617(SSX1 Antibody) at dilution of 1:300