

SSX4 antibody



Catalog Number: 115619

Product name

SSX4 antibody

Immunogen

[Human SSX4 Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

SSX4 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue. Positive WB detected in MCF7 cells, A549 cells, K-562 cells. Observed molecular weight by Western-blot: 22 kDa

Preparation

This antibody was obtained by immunization of SSX4 recombinant protein (Accession Number: XM_017029759). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:115619(SSX4 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human testis using Catalog No:115619(SSX4 antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemical of paraffin-embedded

SSX4 antibody



Catalog Number: 115619

human testis using Catalog No:115619(SSX4 antibody) at dilution of 1:50 (under 10x lens)