

# ST3GAL2 antibody



Catalog Number: 115622

## Product name

ST3GAL2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ST3GAL2 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, HL-60 cells, human colon tissue, human liver tissue, human spleen tissue, K-562 cells. Observed molecular weight by Western-blot: 55kd

## Preparation

This antibody was obtained by immunization of ST3GAL2 recombinant protein (Accession Number: NM\_006927). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:115622(ST3GAL2 antibody) at dilution of 1:500